

SemiFlex™ Distal Extremity Array



Model 165
GE Signa® LX (4 CH)
Model 148
GE Excite® (8 CH)

Description

This 4 to 11 element array (4 of them paired) is dedicated to imaging of the wrist as well as small extremities. This coil conveniently slips onto the wrist or extremity like a bracelet. This is a receive-only phased-array coil.

Suggested Clinical Applications

- Wrist and distal extremity imaging as allowed by the MRI system

Product Specifications

- **Field Strengths:**
1.0T, 1.5T, 3.0T
- **Number of Elements:** 4 to 11
- **Weight:** 2 lbs (0.9 kg)
- **Ship Weight:** 5 lbs (2.3 kg)
- **Cable Length:** 65 in (165 cm)
- **Quality Factor (Q):**
Unloaded: 600 @ 64 MHz
Loaded: 44 @ 64 MHz
- **Physical Size:**

Size	Dia	Length
165	9 cm	18 cm
148	9 cm	18 cm

Warranty/Repair Policy

Product is warranted to be free from defects in materials and workmanship for a period of two (2) years from date of purchase. Typical turn-around time for repair/replacement of defective part(s) is 48 hours after receipt of defective product. Please see the published warranty for details.

System Compatibility

GE Excite®, GE Excite® HD, GE Excite® HDx, GE Signa® LX (4CH), and others. Contact your ScanMed sales representative for details.

Quality Assurance

All products are designed and manufactured under an ISO 13485/9001 registered quality system.



Manufacturer

ScanMed of Resonance Innovations LLC
Authorized Distributor
Radiology OneSource
coils@radiologyonesource.com




+1.615.577.1726
Franklin, TN 37067 USA
www.radiologyonesource.com

SemiFlex™ Distal Extremity Array



Model 165
GE Signa® LX (4 CH)
Model 148
GE Excite® (8 CH)



Manufacturer

ScanMed of Resonance Innovations LLC

Authorized Distributor

Radiology OneSource

coils@radiologyonesource.com

Radiology
OneSource

+1.615.577.1726

Franklin, TN 37067 USA

www.radiologyonesource.com

SemiFlex™ Distal Extremity Array



Model 165
GE Signa® LX (4 CH)
Model 148
GE Excite® (8 CH)



Manufacturer

ScanMed of Resonance Innovations LLC

Authorized Distributor

Radiology OneSource

coils@radiologyonesource.com



+1.615.577.1726

Franklin, TN 37067 USA

www.radiologyonesource.com