



# 1999 GE LIGHTSPEED QXI QUAD



## 1999 GE LightSpeed QXi Quad

Two flat screens  
Octane computer with S DAS  
Generator 140KV at 400MA 56 KW  
Tube last changed August 27<sup>th</sup> 2007 MAS 65,662,668  
Image Capacity 9 GB  
Scan time 0.8  
Reconstruction Time 2 Seconds  
Gantry Count: 409,000,505 Rotations 12,954,802  
512x512 Image Matrix  
MOD and CDRW

Software on Scanner: Software 10.5\_2.8.21\_H1.3M4  
Option, Smart Prep, Connect Pro, Reformat

AW 4.1 WORKSTATION with HP8000 software AW4.1\_04.6\_EXT  
Software options on workstation: Dental plus, Mini viewer, Navigator, Data export 2 web, Add.  
Sub, compare, Fast VR, Reformat, 3D MIP, BW Colon, BWL spine, Cir Aorta 1, Filmier, Mozilla,  
Protocol Manger, Volume rendering.

Ref. 22080

Professionally De-Installed, Crated and Available Immediately

Switzerland

**For more information, please contact us by email at  
[info@radiologyonesource.com](mailto:info@radiologyonesource.com) or call +1.615.577.1717**