



# 2003 GE INFINIA HAWKEYE 4 SPECT/CT

## 2003 GE Infinia Hawkeye 4 SPECT/CT

Camera Infinia II Hawkeye 4

Version 2.104.119.013

Xeleris Version 2.12.20

Dual Monitors

OSI, Recon-jhu, xeleris 2.0, QPS, QGSpect, IOA-6gb, ECTB-ul, CTmrdeconscp, filt3d, rec.base

Camera installation : December 2003

Hawkeye upgrade : April 2008

Tube replacement: : June 2009

Use of the RX tube : 1 or 2 exam/day

Collimators: 1 pair of LEHR

1 pair of collimators especially manufactured by GE for the site ( specifications are intermediate between General Purpose and High Resolution )

1 pinhole collimator

Available August 2010 in Europe

Ref # 24613

Please contact [Debbie van der Touw](#) for more information.

[info@radiologyonesource.com](mailto:info@radiologyonesource.com) or call +1.615.577.1717